

Activity measurement of trypsin using a fluorescence peptide substrate

Introduction

Hydrolysis reaction is caused by making protease act on METHYLCOUMARIN-AMIDE (MCA) of peptide substrate, and 7-AMIDO-4-METHYLCOUMARIN (AMC) isolated.

This isolated AMC, fluorescence becomes the maximum by wavelength 440 nm. Protease activity can be measured using fluorescence spectrophotometer. We introduce the example which performed activity measurement of trypsin using a fluorescence peptide MCA substrate.




Fig. 1 The hydrolysis reaction of MCA by protease

Keywords: Kinetics, Enzyme activity, Lineweaver-Burk, Michaelis-Menten

Measurement system

FP-8300 Spectrofluorometer
 STR-812 Water thermostatted cell holder with stirrer
 CSP-829 Sample compartment lid with syringe port
 MCB-100 Mini Circulation Bath*1)
 VWKN-772 Kinetics Analysis Program

*1) The temperature of a circulation bath is set as 37 degrees by all the measurement.

Samples

Intensity standardization sample: 50 $\mu\text{mol/L}$ AMC solution

Enzyme solution: 10 nmol/L Trypsin bovine pancreas typeVIII, 50 mmol/L Tris-HCl, 0.15 mol/L NaCl, 1.0 mmol/L CaCl_2 , 0.1 mg/mL BSA

Substrate solution: Boc-Gln-Ala-Arg-MCA solution


(The concentration after mixture is adjusted to 0.5, 1, 2.5, 5, 10, 20, 40 $\mu\text{mol/L}$)

Concentration for adjustment. [$\mu\text{mol/L}$]	240	120	60	30	15	6	3
The last concentration. [$\mu\text{mol/L}$]	40	20	10	5	2.5	1	0.5

Measurements

(1) Fluorescence-spectrum measurement of AMC

Excitation and the Fluorescence spectrum of 50 $\mu\text{mol/L}$ AMC were measured (Fig. 2). It turns out that the fluorescence maximum wavelength from this result is 440 nm.



[Measurement parameters]

Ex wavelength: 360 nm

Em wavelength: 440 nm

Ex bandwidth: 5 nm

Em bandwidth: 10 nm

Response: 0.5 sec

Sensitivity: 200 V

Data interval: 1 nm

Scan speed 500 nm/min

Fig. 2 Fluorescence spectrum of AMC

copyright©JASCO Corporation

(2) The vertical axis is changed into concentration from fluorescence intensity

Intensity standardization was performed in order to change the vertical axis into the numerical value equivalent to concentration.

50 $\mu\text{mol/L}$ AMC solution of 0.5 mL is dropped at enzyme solution 2.5 mL. Fluorescence intensity of last concentration of 8.333 $\mu\text{mol/L}$ AMC solution was set to 8333 $\mu\text{mol/L}$ AMC solution.

(3) Enzyme activity measurement

Substrate solution 0.5 mL of each concentration was dropped at enzyme solution 2.5 mL, and time course measurement of the fluorescence intensity of the isolation AMC was performed to it. A result is shown in Fig. 3.




Fig. 3 The pursuit result of an enzyme reaction process

[Measurement parameters]

Ex wavelength: 360 nm	Em wavelength: 440 nm
Ex bandwidth: 5 nm	Ex bandwidth: 10 nm
Data interval: 0.1 sec	Response: 0.1 sec
Scan speed: 200 V	

Analysis

[Kinetics Analysis] In quest of each initial velocity, the Lineweaver-Burk plot was performed using the program from inclination of the time variation data of each substrate concentration (Fig. 4).

$K_m=5.99$ and $V_{max}=35270 \text{ nmol/L}\cdot\text{min}^{-1}$ were obtained from this result.




Fig. 4 Lineweaver-Burk plot