

**Fluorescence Spectrophotometer**

***FP-8000 series***

UV/CD Division  
JASCO Corp



JASCO

# Introduction of New Instruments

JASCO started in 1958 as a company manufacturing  
IR dispersive spectrophotometers

Since then, several other instruments have been  
developed for commercial sales:

UV-Vis(1961)  
CD/ORD(1964)  
Polarimeter(1965)  
Spectrofluorometer(1967)  
Laser Raman(1969)  
HPLC(1972)  
Cell sorter/Flow cytometer(1982)  
SFC(1985)  
Cell fusion(1986)  
Gas permeability tester(1987)  
Ellipsometer(1989)  
NUVOM(1991)  
Ca ion analyzer(1991)  
Hi-SOR-CD(1998)  
SNOM(2000)  
VUV(2000)

The JASCO logo is located in the bottom right corner of the slide. It consists of the word "JASCO" in a stylized, green, sans-serif font.

# History of JASCO Spectrofluorometers




1967 FP-1  
1968 FP-2  
1969 FOM-1  
1970 FP-3  
1970 FL-10  
1971 FOM-2  
1972 FP-4  
1975 FP-100


1976 FP-550  
1976 SFP-3-1  
1980 FP-550A  
1985 FP-770  
1989 FP-777  
1996 FP-750  
1997 FP-715  
2000 FP-6600, 6500,  
6300, 6200, 6100



JASCO

# Diverse Demands of the Market





# Requirements for Life Science Applications

- Micro-sampling
- Automation (Auto-sampler, Micro-plate reader)
- Fluorescence depolarization (anisotropy)
- Auto-titration
- Stopped flow

# Requirements for Advanced Materials Research

- Accurate quantum yield determination
- Accurate spectral correction
- Solid or powder samples
- Micro particle samples
- Luminous materials (White LED & ...)
- Phosphorescence (Organic EL)
- Electric-field inducing Abs. (Organic EL)


# Requirements for NIR Fluorescence Measurements

- Carbon nano-tube analysis
  - Chirality determination of CNT
- NIR & deep red fluorescent dye
  - High sensitivity measurement in crude sample

# Improvement Requests For Next FP Series

- Improvement of noise to signal ratio
- A precise excitation light monitor
- Improvement in user-friendliness
  - Suppression high-order diffraction light
  - Basic and advanced parameters setup
  - Simple and accurate spectrum correction
- Enhanced accessories

## Development Objectives



- 4 models: The successors of FP8200/8300/8500/8600.
  - The fundamental performance of each model's monochromator will be improved.
  - High speed excitation shutter.
  - Auto gain ~ Expand dynamic range.
  - Spectra Manager 2/CFR, Automatic accessories recognition, USB2.0 for data communication.
  - Enhanced accessories.
  - iRM with color LCD
- 

# Instrument Outline

- **Series products by function, performance & applications**
- **Enhancements for performance, S/N, speed, dynamic range, etc.**
- **Improvement in optics and signal processing**
- **Artifact free spectra by higher-order diffraction cut-filters, accurate Excitation monitoring, etc.**
- **Instrument appearance**

# FP-8000 series line

## **FP-8200 High performance routine model**

- Wide dynamic range over 6 orders of magnitude fluorescence emission
- Spectra free from higher order diffraction light (option)

## **FP-8300 Enhanced model for Bio-science**

- Updated application systems for Bio-analysis applications (Microplate reader, Stopped flow, Fluorescence depolarization, Titration)
- Wide dynamic range over 6 orders of magnitude
- Spectra free from higher order diffraction light (standard)
- Expansion to materials analysis (Phosphorescence, Integrating sphere)

## **FP-8500 Evolved for advanced materials evaluation**


- Highest sensitivity
- Highest scan speed
- Wide dynamic range over 6.5 orders of magnitude
- Spectra free from higher order diffraction light (standard)

## **FP-8600 For evaluation of new technology materials**

- Covering NIR region
- Highest scan speed
- Spectra free from higher order diffraction light (standard)

The JASCO logo is located in the bottom right corner of the slide. It consists of the word "JASCO" in a stylized, green, sans-serif font.

## Updated Appearance – Contemporary Design



***FP-8500***

**JASCO**


# Instrument Dimensions





*FP-8200*



*FP-8300*



*FP-8500/8600*



# New Functions and Features

- High S/N
- Wide dynamic range
- High scan speed
- Improved wavelength resolution
- Higher order diffraction cut filter
- Improved NIR performance
- Improved spectral correction
- Enhanced phosphorescence
- Non-contact IQ accessory
- USB communication
- Spectra Manager II & CFR
- iRM-900
- Updated accessories

*FP-8200/8300/8500/8600*

*FP-8200/8300/8500/8600*

*FP-8200/8300/8500/8600*

*FP-8200/8300*

*FP-8200/8300/8500/8600*

*FP-8600*

*FP-8200/8300/8500/8600*

*FP-8300/8500/8600*

*FP-8200/8300/8500/8600*

*FP-8200/8300/8500/8600*

*FP-8200/8300/8500/8600*

*FP-8200/8300*

*FP-8200/8300/8500/8600*

**JASCO**

# Feature 1: Improvement in S/N



The most important feature in Fluorescence

## ***FP-8000***

***FP-8200***            1:380 (SBW=5nm)

***FP-8300***            1:680 (SBW=5nm)

***FP-8500***            1:1200 p-p (SBW=5nm)



***FP-8600***            1:600 (EM SBW=10nm)



**JASCO**

# Feature 1: Improvement in S/N

Water Raman [Peak Signal]/[Peak Noise], Ex 350nm  
**FP-8500: >1200 (RMS), >300 (p-p)**





- Advanced AD enables rapid sampling at the rate 0.2 $\mu$ sec for emission data, and new high speed signal processing by FPGA contributes to considerable improvement of the efficiency for signal acquisition
- P-P S/N is calculated according to JIS K0120

JASCO

# Feature 1: Improvement in S/N

Water Raman [Peak Signal]/[Baseline Noise], Ex 350nm  
**FP-8500: >5000 (RMS)**



- S/N ratio at baseline is much larger than at the Raman peak.
- Noise evaluated at baseline(450nm) is much smaller than at the Raman peak.
- Noise is roughly proportional to signal intensity.

JASCO


## **Feature 2: Expansion of Dynamic Range**

# **New Auto-Gain & Auto-SCS**


**which provides a remarkable  
expansion of fluorescence  
dynamic range**

*JASCO*

## Feature 2: Expansion of Dynamic Range



**A single spectrum collection can provide reliable acquisition of a very small peak together with a coexisting large peak. This feature is very useful for the determination of fluorescence quantum yield.**




JASCO

# Feature 2: Expansion of Dynamic Range

Auto-SCS allows measurement from sub-picomol. to micro-mol, without manually changing the instrument “Sensitivity” setting.

**FP-8200**      **Fluorescein**




**$5 \times 10^{-13} \sim 1.5 \times 10^{-6}$  mol/L**  
**Linearity over 6 orders of magnitude**

$$\log(Y) = A \times \log(X) + B$$

A = 1.007  
B = 9.850  
C.C. = 0.9999  
Std.Err. = 3.469E-008

**FP-8500**      **Fluorescein**



**$5 \times 10^{-13} \sim 1.8 \times 10^{-6}$  mol/L**  
**Linearity over 6.5 orders of magnitude**


$$\log(Y) = A \times \log(X) + B$$

A = 1.023  
B = 9.855  
C.C. = 0.9999  
Std.Err. = 5.585E-008

JASCO

# Feature 2: Expansion of Dynamic Range

Several steps is needed to cover 6 orders in previous model.



## Feature 3: Improvement in Scan Speed

### ***FP-8000***

<b><i>FP-8200</i></b>	20000nm/min
<b><i>FP-8300</i></b>	20000nm/min
<b><i>FP-8500</i></b>	60000nm/min
<b><i>FP-8600</i></b>	Ex:60000nm/min Em:120000nm/min

Advanced sine-bar driving screw and new monochromator motor drive results in a much faster scan speed for all instruments.

JASCO

## Feature 4: Improvement in Wavelength Resolution

***FP-8000***

***FP-8200***                      2.5nm


***FP-8300/8500***                1nm

JASCO


## Feature 5: Higher Order Diffraction Cut Filters

Fluorescence instruments have always suffered from higher order diffraction effects.

Ex spectrum of anthracene





Em Spectrum of Orange plate




JASCO

# Feature 5: Higher Order Diffraction Cut Filters



## Feature 5: Higher Order Diffraction Cut Filter



**The new 'cut filter' capability removes peaks originating from higher order diffracted light to provide simple, reliable spectral acquisition and analysis. This system works quite well for 3D fluorescence measurements to reveal peaks usually hidden by the high order diffraction peak. The cut filter method also works quite well for quantum yield determination by providing measurements of accurate spectra.**

JASCO


## Feature 6: Improvements for Red Wavelengths

Rare-earth compounds, such as Eu, are often used in luminous materials such as LEDs which are evaluated by fluorescence in the deep red wavelengths.

To satisfy needs for these applications, the FP-8500/8600 instruments are equipped with a focus adjustment mechanism to maintain high wavelength resolution over a wide spectral range.

## Feature 6: Improvements for Red Wavelengths


### Em spectrum of a ruby



JASCO

## Feature 7: High Accuracy Spectrum Correction

The standard fluorescence spectral correction method was established long ago and described in many textbooks by using a Rhodamine B standard and a calibrated W light source. It is generally accepted that this spectral correction method may not provide the most accurate spectral correction, but it is abided by the inevitable that fluorescence spectrophotometers has such characteristics by its nature.




The evaluation of luminous compounds and determination of quantum yield for advanced materials research requires a more accurate spectral correction method. Jasco was determined to provide a more accurate spectral intensity correction method and the tools required for fluorescence correction.

## Feature 7: High Accuracy Spectral Correction

To provide greater accuracy, Jasco developed a calibrated light source – D2 and a calibrated detector for the FP-8000 series, which will cover the full wavelength range with directly observed correction data.

FP-8000 Correction Method



Rhodamine B



Calibrated light source - D2



SID-844 Calibrated detector



JASCO

## Feature 8: Enhancements for Phosphorescence

Organic Electro-Luminescence (EL) is a new application for phosphorescent materials.

To analyze this specific type of sample, we developed a rotating chopper mechanism to enable measurements of phosphorescent materials with a lifetime shorter than 1 msec.


- Rotation frequency can be selected from 1200, 600, 300, 150 or 75 rpm
- Rotating chopper allows quick phosphorescence spectra measurement
- Expand phosphorescence measurements to many modes:  
([Spectra], [Calibration]/[Analysis], [Fixed wavelength], [Time course], [Phosphorescence lifetime])



JASCO


## Feature 8: Enhancements for Phosphorescence

### Phosphorescence lifetime of Benzophenone




Range : 27.4 - 43.3 msec  
Life time : 5.40 msec  
Std. err. : 0.001799  
Regression:  $Y = -0.185028 \times X + 11.9925$   
C.C. : -0.9996

### Fluorescence of benzene





### Phosphorescence of benzene



JASCO

## Feature 9: Non-contact IQ Accessories

Utilizing a non-contact type RFID tag




JASCO

# Feature 10: Other features



Spectra Manager Version2,  
Spectra Manager CFR  
(FP-8200/8300/8500/8600)




iRM-900  
(FP-8200/8300)




USB connection





# Introduction of Individual Models



# ***FP-8500***

**Highest performance in the world  
Enhanced for solid sample &  
advanced materials analysis**



**JASCO**

## ***FP-8500***



- **Highest sensitivity**  
(>1200 RMS peak, >5000 RMS base)
- **Highest scan speed**  
(Max: 60,000 nm/min)
- **Dynamic range of 6.5 orders of magnitude**
- **Reduction of higher order diffraction peaks**


Evaluation of advanced materials: white LEDs, organic EL, etc., increased the requirements for fluorescence instruments.

To meet these applications, we have enhanced the spectral correction tools, integrating sphere, and other accessories together with basic functions and performance; the S/N, scan speed, dynamic range, and other enhancements.

**JASCO**

# FP-8500

## Wide dynamic range



Linearity over 6.5 orders of magnitude  
 $5 \times 10^{-13} \sim 1.8 \times 10^{-6}$  mol/L of fluorescein

$$\log(Y) = A \times \log(X) + B$$


$$A = 1.023$$

$$B = 9.855$$

$$\text{C.C.} = 0.9999$$

$$\text{Std.Err.} = 5.585 \text{ E-008}$$


## Highest sensitivity




Peak to Peak	352
RMS (peak)	1534
RMS (base)	8463

## Highest scan speed

### Quinine sulfate





# ***FP-8600***

To meet new applications, such as carbon nano-tubes, NIR fluorescent dyes, up-conversion of fluorescence glasses, and other NIR fluorescent materials to be developed.



JASCO

## ***FP-8600***



**Carbon nano-tubes, NIR fluorescent dyes, up-conversion of fluorescent glasses.**


**New materials such as these will require NIR analysis.**

**We have enhanced the performance and functions to meet these applications.**


- **NIR model for evaluation of new materials  
(Ex200nm~850nm, Em200nm~1010nm)**
- **Highest scan speed  
(Max: Ex 60,000 nm/min, Em 120,000 nm/min)**
- **Reduction of higher order diffraction peaks**

**JASCO**


**FP-8600**



## Ruby



**JASCO**



# ***FP-8300***

Targeted for Bio- and Life Science  
applications



JASCO

## ***FP-8300***




- **Enhanced for Bio-applications**
- **Wide dynamic range over 6 orders of magnitude**
- **Reduction of higher order diffraction**
- **Expanded to solid samples (Phosphorescence, integrating sphere)**
- **High sensitivity (>680 RMS peak, >2800 RMS base)**
- **High speed (Max: 20,000 nm/min)**

**Bio-analysis is the most common application for fluorescence analysis. We have drastically improved the performance for the FP-8300, and developed enhanced accessories dedicated for these applications, such as the micro-plate reader, stopped flow, fluorescence anisotropy and auto-titration.**


**JASCO**

# FP-8300


## Fluorescence depolarization anisotropy




## Auto-titration




## Stopped flow



## pH Phosphorescence



JASCO



# ***FP-8200***

Routine fluorescence analysis



JASCO

## ***FP-8200***



We have equipped the FP-8200, our routine model, with valuable capabilities, such as Auto Gain, Auto SCS and an optional cut-off filter assembly, for easy, quick, and reliable analyses. We have also improved the sensitivity, scan speed and resolution, among other items.


- Wide dynamic range over 6 orders of magnitude
- Reduction of higher order diffraction peaks (option)
- High sensitivity (>380 RMS peak, >1600 RMS base)
- High speed (Max 20,000 nm/min)

**JASCO**

# FP-8200

Wide dynamic range

Fluorescein



Linearity over 6 orders of magnitude:

$5 \times 10^{-13} \sim 1.5 \times 10^{-6} \text{ mol/L}$

$$\log(Y) = A \times \log(X) + B$$

$$A = 1.007$$


$$B = 9.850$$

$$\text{C.C.} = 0.99995$$

$$\text{Std.Err.} = 3.469\text{E-}008$$

High sensitivity

Diluted fluorescein samples



- $1 \times 10^{-10} \text{ mol/L}$
- $5 \times 10^{-11} \text{ mol/L}$
- $1 \times 10^{-11} \text{ mol/L}$
- $5 \times 10^{-12} \text{ mol/L}$
- $1 \times 10^{-12} \text{ mol/L}$
- $5 \times 10^{-13} \text{ mol/L}$

JASCO

# Specifications



Model	FP-8200	FP-8300	FP-8500	FP-8600	
Light source	150W Xe lamp with shielded lamp housing				
Photometric system	Photometric ratio system using monochromatic light to monitor the intensity output of Xe source				
Cut-off filter assembly for higher order diffraction	Option	Standard			
Sensitivity* <sup>1</sup> (RMS)	Peak* <sup>2</sup>	380:1	680:1	1200:1	600:1
	Base* <sup>3</sup>	1600:1	2800:1	5000:1	2500:1
Wavelength range	Ex: Zero order , 200-750nm (Optionally -850 nm) Em: Zero order , 200-750nm (Optionally -900 nm)	Zero order , 200-750nm (Optionally -900 nm) on both Ex and Em	Zero order , 200-750nm (Optionally -850 nm) on both Ex and Em	Ex: Zero order , 200-850nm Em: Zero order , 200-1010nm	
Spectral bandwidth	2.5-20nm on both Ex and Em	1-20nm on both Ex and Em	1-20, L5, L10nm on both Ex and Em	Ex: 1-20, L5, L10nm Em: 2-40, L10, L20nm	
Wavelength scan speed	20-20000nm/min		10-60000nm/min	Ex: 10-60000nm/min Em: 20-120000nm/min	
Resolution	2.5nm	1.0nm		Ex: 1.0nm Em: 2.0nm	
Wavelength accuracy	±2.0nm	±1.5nm	±1.0nm	Ex: ±1.0nm Em: ±2.00nm	
Gain	High, Medium, Low, Very Low, Auto				
Auto Gain	Standard				
Start button	Standard				
Dimensions	490(W) × 545(D) × 270(H)mm	520(W) × 545(D) × 270(H)mm	570(W) × 545(D) × 270(H)mm		
Net weight	33.6kg	36kg	39kg		

\*<sup>1</sup>: Minimum Signal to Noise ratio of Raman band of water, excitation 350nm, band width Ex 5nm Em 5nm (FP-8600: Ex 5nm Em 10nm), response 2s

\*<sup>2</sup>: Noise is measured on the Raman peak

\*<sup>3</sup>: Noise is measured on the baseline



JASCO

# *Accessories*

Important instrument interfaces for a variety of sample analyses requirements



JASCO

## Accessories for various volumes and shapes of samples

Model:	Name:	Available for:
FMH-801	3 mm micro cell jacket	FP-8200 / 8300 / 8500 / 8600
FMH-802	5 mm micro cell jacket	FP-8200 / 8300 / 8500 / 8600
FUV-803	Abs. measurement cell block	FP-8200 / 8300 / 8500 / 8600
FHM-804	High sensitivity cell block	FP-8200 / 8300 / 8500 / 8600
FSA-805	Cell block for triangular cell	FP-8200 / 8300 / 8500 / 8600
FSA-806	Cell block for rectangular cell	FP-8200 / 8300 / 8500 / 8600
FDA-808	Solid sample holding block	FP-8200 / 8300 / 8500 / 8600
FLH-809	Film holding block	FP-8200 / 8300 / 8500 / 8600
FPA-810	Powder sample cell block	FP-8200 / 8300 / 8500 / 8600



FMH-801/802



FUV-803



FHM-804



FSA-805



FSA-806



FDA-808



FLH-809



FPA-810

JASCO

## Accessories for various volumes and shapes of samples (2)

Model:

SAF-850

SAF-851

CTS-855

OBF-832

EFA-833

Name:



One-drop measurement unit



One-drop measurement unit



Coumarin Measurement Unit


Optical fiber unit

Epi-fluorescence unit

Available for:

FP-8200

FP-8300 / 8500 / 8600



SAF-850



SAF-851




OBF-832



EFA-833

JASCO


# EFA-833 Epi-fluorescence Accessory



- L1, L2, L3: **Lenses 1, 2, 3**
- M1, M2, M3: **Mirrors: 1, 2, 3**
- A: **Aperture**
- ES: **Quartz window**
- L4, L5: **Lenses 4, 5**
- M4, M5: **Mirrors 4, 5**

Optics with mirrors and lenses improves Ex. light energy markedly, especially in far UV region.

Ex. intensity observed with Conc. RhodamineB



- EFA-833 Epi-fluorescence unit
- FDA-808 Solid sample block

The EFA-833 provides approximately 20% greater Ex. intensity than the FDA-808.

JASCO

# Accessories for Temperature Control

Model:	Name:	Available for:
CTH-807	Water thermostatted cell block	FP-8200 / 8300 / 8500 / 8600
STR-811	Water thermostatted cell holder with stirrer	FP-8200
STR-812	Water thermostatted cell holder with stirrer	FP-8300 / 8500 / 8600
FCT-816	 Water thermostatted 4-position cell changer	FP-8200
FCT-817	 Water thermostatted 8-position cell changer	FP-8300 / 8500 / 8600



CTH-807



STR-811



STR-812



FCT-816



FCT-817

# Accessories for Temperature Control

Model:	Name:	Available for:
EHC-813	Peltier thermostatted single cell holder (air-cooled)	FP-8300 / 8500 / 8600
ETC-814	Peltier thermostatted single cell holder (water-cooled)	FP-8200
ETC-815	Peltier thermostatted single cell holder (water-cooled)	FP-8300 / 8500 / 8600
PCT-818	Automatic 4-Position Peltier cell changer (water-cooled)	 FP-8300 / 8500 / 8600
PMU-830	Liq. N <sub>2</sub> cooling unit	FP-8300 / 8500 / 8600
CSH-831	 Cryostat holder	FP-8300 / 8500 / 8600
HPC-836	 High temperature powder cell unit	FP-8300 / 8500 / 8600



EHC-813



ETC-814



ETC-815



PCT-818



PMU-830



Temp. controller




CSH-831



HPC-836

- PCT-818 is adaptable to FP-8300
- CSH-831 is included in series
- HPC-836 is controllable from PC

# Accessories: Automatic Polarizer




Model:  
**FDP-837**



Name:  
**Automatic polarizer unit**

Available for:  
**FP-8300 / 8500 / 8600**



**Driving motor**



**Top: Quartz window  
Bottom: Polarizer**

**We have separated the polarizer assembly from the system, and enhanced its flexibility, which allows the use of a thermostatted cell holder with a circulating bath instead of a dedicated Peltier cell holder. The new polarizer system is also for the FP-8300. This new accessory can also be applied to the fluorescence polarization measurement of solid samples.**



# Accessories ~ Auto-sampler and Sipper

Model:	Name:	Available for:
SHP-819	Peristaltic sipper	FP-8200
SHP-820	Peristaltic sipper	FP-8300 / 8500 / 8600
<b>QFS-821</b>	<b>Vacuum sipper</b>	<b>FP-8200</b>
QFS-822	Vacuum sipper	FP-8300 / 8500 / 8600
FSC-823	Micro flow cell holder	FP-8200
FSC-824	Micro flow cell holder	FP-8300 / 8500 / 8600
<b>FMP-825</b>	<b>Automatic microplate reader</b>	<b>FP-8300 / 8500 / 8600</b>
ASU-800	Auto sampler unit	FP-8200 / 8300 / 8500 / 8600
ASP-849	Syringe pump unit	FP-8200 / 8300 / 8500 / 8600



SHP-819



QFS-821



FSC-823



FMP-825



ASU-800



ASP-849



SHP-820



QFS-822



FSC-824

- New vacuum sipper was developed for FP-8200
- FMP-825 was modified



## FMP-825 Automated Microplate Reader

Faster measurements are obtained by modifying the method for positioning of the micro-well plate and data acquisition.  
We have also modified the temperature controller, and added as optional functions.




JASCO

# Accessories for Quantum Yield Determination

Model:	Name:	Available for:
ISF-834	60mm $\Phi$ Integrating sphere unit	 <b>FP-8300</b> / 8500 / 8600
ILF-835	100mm $\Phi$ Integrating sphere unit	 <b>FP-8300</b> / 8500 / 8600
ILFC-847	Cooled 100mm $\Phi$ Integrating sphere unit	 <b>FP-8300</b> / 8500 / 8600



ISF-834



ILF-835




ILFC-847

•Adaptable to FP-8300

JASCO

# Accessories for Spectral Correction

Model:	Name:	Available for:
ESC-842	 Calibrated light source – WI	FP-8200 / 8300 / 8500 / 8600
ESC-843	 Calibrated light source – D2	FP-8200 / 8300 / 8500 / 8600
SID-844	 Calibrated detector	FP-8200 / 8300 / 8500 / 8600



ESC-842



ESC-843



SID-844

- **ESC-842 is now available for the FP-8200 as corrected spectra are indispensable for publishing of research papers.**
- **The calibrated light source units were modified to fit the std. cell holder so the sources can be used directly for proper operation**
- **ESC-843 and SID-844 are newly developed spectral correction tools.**

# Accessories for Kinetics and Reaction Analyses

Model:

SFS-852 / 853/ 854

ATS-826

ATS-827

CSP-828

CSP-829

Name:



Stopped flow system



Automatic titration unit



Automatic titration unit

Lid with syringe port

Lid with syringe port

Available for:

FP-8300 / 8500 / 8600

FP-8200

FP-8300 / 8500 / 8600

FP-8200

FP-8300 / 8500 /8600



SFS-852/853/854



Cell



ATS-826



ATS-827



CSP-828



CSP-829

# Validation Accessories and Depolarizer

Model:

VDK-840

VDK-841

FSP-838

Name:

 Validation kit 1

 Validation kit 2

 Depolarization plate

Available for:

FP-8200 / 8300

FP-8200 / 8300 / 8500 / 8600

FP-8200 / 8300 / 8500 / 8600

## VDK-840

Testing kit for stray light and spectra correction.  
(FP-8500/8600 included as standard)

**Constitution:**  
Rhodamine B  
(bottled, with pipet)  
Triangle cell w/  
PTFE stopper

## VDK-841


Filter set (with certificate) for testing stray light

**Constitution:**  
WG-305  
L39  
R60



FSP-838

JASCO




# **System control/ Data processing**



JASCO

# Spectra Manager Version 2/CFR

**SM Ver2/CFR  
with  
numerous enhancements**



**JASCO**


# Standard Programs (DS measurement)

Name	Available for:	CFR
Spectrum measurement	2 / 3 / 5 / 6	Yes
Time course measurement	2 / 3 / 5 / 6	Yes
Quantitative calibration	2 / 3 / 5 / 6	Yes
Quantitative measurement	2 / 3 / 5 / 6	Yes
Fixed wavelength measurement	2 / 3 / 5 / 6	Yes
 Abs measurement	2 / 3 / 5 / 6	Yes
Intensity monitor	2 / 3 / 5 / 6	Yes
Validation	2 / 3 / 5 / 6	Yes
Parallel time course	2 / 3 / 5 / 6	Yes
Phosphorescence lifetime	3 / 5 / 6	Yes
3D spectra measurement	2 / 3 / 5 / 6	Yes
Interval scan measurement	2 / 3 / 5 / 6	Yes
Spectral correction	2 / 3 / 5 / 6	Yes


Std. programs (DS measurement) ~ Phosphorescence mode

**These programs are provided with a 'phosphorescence mode'.**

Name	Phosphorescence available for
Spectrum measurement	3 / 5 / 6
Time course	3 / 5 / 6
Quantitative calibration	3 / 5 / 6
Quantitative measurement	3 / 5 / 6
Fixed wavelength	3 / 5 / 6




**Phosphorescence time course measurement can be used for the evaluation of light-accumulating (phosphorescent) pigment, while phosphorescent assays are a more popular application for bio-science.**




JASCO

# Std. programs (DS measurement) ~ Phosphorescence mode



## Std. programs (DS measurement) ~ Time course



### FP-8200

- Data pitch: 0.01 sec ~ 60 sec
- Max measure time: approx. 7 days (600,000 sec/10,000 min/167 hr)

### FP-8300/8500/8600

- [Fluorescence mode] & [Phosphorescence mode], are both available
- [Fluorescence mode]
- Data pitch: 0.05 msec ~ 60 min
- Unit of time: msec, sec, min
- Max measure time: 200 days (288,000 min/4,800hr)


Light-accumulation (phosphorescent) pigments have a life time of minutes to hours. Lifetime can be determined by observing the emission time course data after shutting off the Ex. light.




# Std. programs (DS measurement) ~ 3D spectra




## **New** Simpler parameter settings




## 3D spectra



## Contour map



## Color contour




JASCO


## Std. programs (DS measurement) ~ Abs measurement

Many textbooks say “Fluorescence from samples may provide an artifact to the observed absorbance spectra” without any method for ‘correction’ of this interference.

The FP-8000 series offers a perfect absorbance measuring system: using the dual monochrometers, the higher order diffraction cut filter units for both incident and transmitted light and the use of the synchronized scanning mode with 0 nm wavelength difference. It is now possible to observe accurate absorbance characteristics of organic EL compounds and other fluorescent samples.



Abs. spectrum of Eosin B



JASCO

## Std. programs (DS measurement) ~ Spectral correction

- Calibrated D<sub>2</sub> light source and a calibrated detector unit are new developments. They, together with Rhodamine B and the calibrated WI light source provide the ability to obtain correction data for the full wavelength range for both Ex. and Em. directly.
- These calibrated components are applicable to not only the std. cell holder but also the ISF-834 60mm  $\Phi$  and ILF-835 100mm  $\Phi$  I.S. units.
- The correction program has been modified to be able to operate Ex. and Em. separately.

Ex \ Em	1	2.5	5	10	20
1	○	○	○	○	○
2.5	○	○	○	○	○
5	○	○	○	○	○
10	○	○	○	○	○
20	○	○	○	○	○

Previously, 50 spectra were required to cover all possible parameter settings.



Ex	
1	○
2.5	○
5	○
10	○
20	○
L5	○
L10	○


Em	
1	○
2.5	○
5	○
10	○
20	○
L5	○
L10	○

The correction system has been improved to reduce the requirements to 14 spectra.

JASCO

# Std. programs (DS measurement) ~ Peak search


**Peak search**  provides a rapid method to determine the Ex. and Em. peaks for an unknown fluorophore. This function has been modified to be accessible from each software application, instead of an independent function. This enhances the ease of operation.



## Std. programs (DS measurement) ~ 'Single-point calibration'

The FP-8000 series is equipped with a new function which can calibrate the fluorescence intensity with the observed intensity of a single standard. This function eliminates the need to generate a working calibration curve at the beginning of every fluorescence assay.

For example, the regulation for the “testing method of chlorophylls” as an authorized testing method for water requires a calibration of the instrument by adjusting the fluorescence intensity of 1 mg/mL of fluorescein to 100.



## Std. programs (DS analysis)

Name	Available for:	CFR
Spectrum analysis	2 / 3 / 5 / 6	Yes
Interval data analysis	2 / 3 / 5 / 6	Yes
JASCO Canvas	2 / 3 / 5 / 6	Yes


### Programs in Spectrum analysis:

Name	Available for:	CFR
Phosphorescence lifetime	3 / 5 / 6	Yes
Spectral correction	2 / 3 / 5 / 6	Yes
 <b>Relative quantum yield</b>	2 / 3 / 5 / 6	No
Enzymatic reaction rate	2 / 3 / 5 / 6	Yes
 <b>X unit conversion</b>	2 / 3 / 5 / 6	Yes
 <b>Y unit conversion</b>	2 / 3 / 5 / 6	Yes
 <b>Time axis offset</b>	2 / 3 / 5 / 6	Yes
 <b>Peak normalization</b>	2 / 3 / 5 / 6	Yes

JASCO

## Std. programs (DS analysis) ~ Relative quantum yield


Quantum yield is proportional to the ratio of total emission when both fluorophore have the same total absorbance. A new program has been developed to determine the quantum yield of unknown samples by direct comparison with a standard sample having a known quantum yield.



JASCO

# Std. programs (DS analysis) ~X unit conversion

A new X unit conversion program has been added which does not distort the peak area. This conversion method is physically important when the ordinate has a dimension such as number of photons, energy, etc.




## Std. programs (DS analysis) ~Y unit conversion

A fluorescence photometer provides a spectrum with an ordinate value with the dimension of number of photons. While suited for quantum yield determination, the dimension of 'energy' is needed for color analysis of fluorescent emission. This new function provides such an ordinate conversion without changes to the dimensionless value of 'number of photons'.


For example, one photon has the energy of  $e=h\nu$ , so the sum of  $N$  photons is  $E=N\cdot h\nu=N\cdot hc/\lambda$ .

This function converts the ordinate by multiplying/dividing the factor of wavelength such that  $1\ \mu\text{m}$  is set to be 1.



## Std. programs (DS analysis) ~ Peak normalization


**This function normalizes the scales of multiple spectra for simplified shape comparison. This function is very useful for direct comparison of fluorescent dyes.**



JASCO

## Std. programs (DS analysis) ~ Time axis offset

**This function provides an adjustment for setting of a proper “offset” to the time axis of a time course observation. It is valuable for kinetics analysis to compensate for the time lag from “measurement start” to “reaction start”.**



JASCO

## Std. programs (Driver)



Name	Available for:
FP Control Driver	2 / 3 / 5 / 6
PAC-743 Temperature Controller Driver	2 / 3 / 5 / 6
Temperature Controller(ETC, EHC, PSC)	2 / 3 / 5 / 6
Driver for High Temperature Controller	3 / 5 / 6
Driver for AUS-800 Autosampler	2 / 3 / 5 / 6
Driver for AUS-605 Autosampler	2 / 3 / 5 / 6
Driver for Microplate Reader	3 / 5 / 6
Control driver for the rapid stopped-flow system	3 / 5 / 6
Driver for Automatic Titration Unit	2 / 3 / 5 / 6




# Optional programs ~ DS

Model	Name	Components	Available in	CFR
FWTS-872	Temperature interval	Temperature interval scan	2 / 3 / 5 / 6	Yes
FWTC-873	Dual-wavelength time course	Dual-wavelength time course Intracellular ionic concentration	2 / 3 / 5 / 6	Yes
FWTP-874	Temperature control – Melting	Temperature control/ramping DNA melting analysis Melting temperature calculation	2 / 3 / 5 / 6	Yes
VWKN-775	Kinetics analysis	Kinetics analysis	2 / 3 / 5 / 6	Yes
FWAP-875	Auto-depolarization	Auto-depolarization fixed W.L. Auto-depolarization time course	3 / 5 / 6	Yes
		Auto-depolarization temp. ctrl.		
		<b>Auto-depolarization titration</b>		
		<b>Depolarization analysis</b>		
FWSF-877	Stopped flow	Stopped-flow measurement	3 / 5 / 6	No
		Reaction rate calculation		
		Data accumulation		
FWAT-876	Auto-titration	Automatic titration measurement	2 / 3 / 5 / 6	Yes
FWFC-878	Fluorescent object color	Fluorescent object color	3 / 5 / 6	No
FWLU-879	Luminous color	Luminous color meas./anal.	3 / 5 / 6	No
VWMC-883	Macro command	Macro command	2 / 3 / 5 / 6	No
FWQE-880	Quantum yield calculation	Quantum yield calculation	3 / 5 / 6	No

# *IRM-900*

Modernized color LCD panel



JASCO

# Std. programs (iRM)



## Measurement


## Analysis

Name	Available in	Name	Available in
Spectra	2 / 3	 spectra analysis	2 / 3
Time course	2 / 3	 3D spectra analysis	2 / 3
Quantitative meas.	2 / 3		
Fixed wavelength	2 / 3		
 Abs measurement	2 / 3		
Intensity monitor	2 / 3		
Peak search	2 / 3		
Validation	2 / 3		
 3D spectra	2 / 3		
 Spectral correction	2 / 3		




# Std. programs (iRM) ~ for easier operation

- Analysis separated from measurement with enhanced functions
- Two modes of parameter settings: basic and advanced modes



**Spectrum measurement**



**Spectrum analysis**

## Data processing:

- Baseline correction
- Smoothing
- Peak cut
- Deconvolution
- FFT filtering
- Interpolation, Cut
- Arithmetics
- Differential
- K-K (Kramers-Kronig)
- Peak picking
- Peak height/ratio
- Peak area/ratio
- Peak half-height width
- Spectra difference
- X/Y unit conversion
- Comments editor
- Enzymatic reaction rate
- Spectra

JASCO

# Std. programs (iRM) ~ 3D spectrum


## 3D Spectra measurement

- Up to 50 spectra can be acquired as interval data (jwb)


## 3D spectra analysis

- Display: Color contour / Contour / 3D / Color 3D
- Data processing: Smoothing / Baseline correction / Comment Editing
- Cross section, Peak search


Color contour




Contour



3D



Color 3D




JASCO

# Optional programs ~ iRM



Model	Name	Component	Available in
 FRTC-884	Dual-wavelength time course	Dual-wavelength time course	2 / 3
 F RTP-882	Temperature controlling • Melting	Temperature controlling DNA Melting analysis	2 / 3 2 / 3
FRKN-881	Kinetics	Parallel time course Kinetics analysis	2 / 3
 FRAP-883	Auto-depolarization	Auto-depolarization fixed wavelength Auto-depolarization time course Auto-depolarization temperature controlling Depolarization analysis	3
FRMC-885	Macro command	Macro command	2 / 3





Fluorescence Spectrophotometer

***FP-8000 series***

Application system



UV/CD Division  
JASCO Corp.

JASCO

# Application system

- Optimal combination with accessory and program dedicated to specific application

# Application systems

- Phosphorescence
- Fluorescence depolarization
- Fluorescence color
- Fluorescence quantum yield
- Phosphorescence quantum yield
- Auto measurement
- Stopped flow
- Auto titration
- One drop fluorescence
- Single particle fluorescence

# Phosphorescence


Organic EL is a growing application for phosphorescence. For R&D research, a material is analyzed for a triplet emission, and then modifications are attempted to improve the quantum yield. By obtaining a more efficient excitation to the triplet state, energy losses are minimized and higher efficiency is obtained. Phosphorescence assays are also increasingly popular for bioscience applications.

## System:

- FP-8500 (8300, 8600)
- PMU-830 Liq. N<sub>2</sub> cooling unit



## Lifetime of benzophenone



Range: 27.4 - 43.3 msec  
Life time: 5.40 msec  
Std.error: 0.00180  
Regression:  $Y = -0.1850 \times X + 11.99$   
C.C.: -0.9996

JASCO


# Auto depolarization measurements

Anisotropy is a well established method to study the mobility of a molecule, the molecular size, quality of binding to receptor, etc, by observing the depolarization of fluorescence from molecules excited by linear polarized light.

The polarizer drive unit is now a separate unit and can be combined with any cell holder of the customer's requirement.

Example system:


- FP-8300 (8500)
- ETC-815 Peltier thermostatted cell holder
- FDP-837 Automatic polarizer
- FWAP-875 Auto-depolarization program



Polarizer motor drive unit



Top: Quartz window  
Bottom: Polarizer



Fluorescent dye (DPH) was added to a liposome sample, then the fluorescence polarization was observed by changing the temperature. This experiment reveals a phase transition at approximately 40° C.

JASCO

# Color evaluation of fluorescent materials




Fluorescent materials are increasingly used for many objects, such as safety signs for superior visibility. ASTM established a testing method for luminescent color by using a fluorescence spectrophotometer. The standard method requires the measurement of the 3D spectra, calibration of the Ex. intensity to the standard D<sub>65</sub> source, and accumulation of the Em. intensity for each wavelength for an intensity distribution. The object color calculation is then performed to obtain the chromaticity coordinates for the sample.

## System

- FP-8500 (8300)
- ISF-834 60mm diam. I.S. accessory
- FWFC-878 Fluorescent object color program



Spectra of fluorescent sample under D<sub>65</sub>







JASCO

# Fluorescence quantum yield determination

Fluorescent systems has been changed by the requirement from the evaluation of luminous materials used for white LED, etc. The FP-8000 series instruments provide numerous features for optimized analysis of these materials.

## System

- FP-8500 (8300)
- ISF-834 60mm diam. I.S. acc'y or ILF-835 100mm diam. I.S. acc'y
- FWQE-880 Quantum yield program



$$\Phi = \frac{\text{\# Emitted photons}}{\text{\# Absorbed photons}} = \frac{S_2}{S_0 - S_1}$$



JASCO

# Phosphorescence quantum yield determination

An optimized accessory for phosphorescence quantum yield determination has been designed which includes a dewar with an integrating sphere.

## System

- FP-8500
- IFLC-847 Liq. N<sub>2</sub> cooled 100mm diam. I.S. acc'y
- FWQE-880 Quantum yield program



$$\Phi = \frac{819.276}{4954.31 - 4073.96} = 0.93$$




(Published spec.  $\Phi > 0.8$ )

Exceeds the published specification by one figure in the precision.

JASCO

# Auto-sampler measurements


The standard measurement programs: [Spectrum Measurement], [Quantitative Analysis] and [Fixed Wavelength Measurement] for the PC or iRM includes an dedicated interface for automated measurements using an analysis sequence which can be modified by the user. This software function is usable for all three auto-systems: peristaltic sipper, vacuum sipper, and flow cell holder + syringe pump.




# Stopped-flow experiments

Stopped flow experiments are widely applied to the study of protein denaturation, enzyme reaction, and others. The combination of a rapid and sensitive fluorescence spectrophotometer and a stopped flow unit provides an optimized system.

## Mixing Efficiency




Calc. S.D.  
of ratio



entrance

position


exit



Entrance  
of the cell

position

Exit  
of the cell



Denaturation of Cytochrome C (0.5mg/mL) by 0.1N sulfuric acid. Mixing ratio of 1:1.

JASCO

# Auto-titrator experiments

CD is widely used for protein denaturation studies using pH change or other denaturation compounds. Fluorescence can also be applied for these experiments. JASCO offers a new integrated system incorporating an updated auto-titration unit and a dedicated program which will surpass competitive systems.


## System

- FP-8300 (8200, 8500)
- STR-812 Water Thermostatted cell holder with stirrer
- ATS-827 Automatic titration unit
- FWAT-876 Automatic titration program



pH dependency of fluorescence intensity of BCECF

BCECF Fw:536.5



pH sensitive fluorescent dye


Ex  $\lambda$  : 505 nm

Em  $\lambda$  : 530 nm

Isosbestic point: 442 nm

pKa=6.98

Reagent: 0.005 N NaOH



JASCO

# One drop measurement


“The smallest amount of samples possible” is an emerging requirement, especially for bio-science applications. The success of the “one-drop measurement” for the V-600 series has been extended to the FP-8000 series.

## System

- FP-8200 (8300,8500,8600)
- SAF-850 One-drop meas. acc’y (SAF-851)





SAF-850



## DNA labeled with PicoGreen

A  $\lambda$  DNA sample was labeled with the fluorescent dye PicoGreen, which is a highly selective and sensitive label for the double strand DNA. The one drop experiment demonstrates more than 1000X higher sensitivity than the UV absorbance method.




# Single particle fluorescence

Several researchers have requested the ability to determine the fluorescence quantum yield of a single luminous particle of a sample ranging from 10-100  $\mu\text{m}$  in size. The KBr plate method can provide the ability to perform such an experiment.


## System

- FP-8500
- ILF-835 100mm diam. I.S. acc'y
- KBr plate holder (custom)
- Optional KBr plate sampling tool




Patent pending


Fluorescence spectrum of one grain sample



Effectively usable for phosphorescence measurement



JASCO



**Thank you for attention!**



*JASCO*