

Analysis of Organic Acid in Alcoholic Beverages

Introduction

Organic acids are compounds that produce tastes such as sourness and deliciousness. They also have some significant features such as oxidation resistant and antibacterial activity. To detect organic acids by HPLC, some methods have been prepared, such as detection with low wavelength of UV, direct detection using Refractive Index detector or electric conductivity detector, and detection with using a pH indicator. The direct detection method, however, has difficulties to separate the target and foreign substances so that the quantification of target substances is difficult. High mobile phase background levels and low selectivity are also problematic. While pH indicator detection can be used, higher selectivity and sensitivity analysis can be achieved through visible absorbance detection by a post column method.

In this application, the analysis of organic acids in alcohol is implemented using a post column method and BTB as pH indicator.


The detection wavelength of UV detector is 445nm.

Keyword: Organic acids, alcohol, pH indicator, BTB, UV-Vis detector, Shodex RSpak KC-811

Experimental Condition

Column:	Shodex RSpak KC-G + KC-811 x 2 (6.0 mmI.D. x 50 mmL) (8.0 mmI.D. x 300 mmL x 2)	Eluent:	2mM Perchloric acid
		Reagent:	0.2mM BTB + 15 mM Disodium hydrogen phosphate
Eluent flow rate:	1.0 mL/min	Column temp.:	60°C
Reagent flow rate:	1.5 mL/min	Injection volume:	50 µL
Wavelength:	445 nm	Sample:	Alcoholic beverages
Standard:	14 organic acids, 1 mmol/L each		

Schematic diagram



- 1: Eluent
- 2: Degassing unit
- 3: Eluent pump
- 4: Autosampler
- 5: Gard column (RSpak KC-G)
- 6, 7: Column (RSpak KC-811)
- 8: Column oven
- 9: Reagent
- 10: Reagent pump
- 11: Organic acid reaction unit
- 12: UV/VIS detector
- 13: Chromatography data system (ChromNAV Ver.2)

Pretreatment

Rice wine: Rice wine is filtered by 0.45 µm membrane filter

Red wine: After the red wine is diluted ten times by ultrapure water, the solution is filtered by 0.45 µm membrane filter

Plum wine: After the plum wine is diluted ten times by ultrapure water, the solution is filtered by 0.45 µm membrane filter

Results

Figure 1 shows chromatograms of organic acids mixture solution. This alcohol sample contains not only twelve organic acids mixture standard sample but also phosphoric acid and tartaric acid. Therefore, the standards of the two substances are also measured.




Fig.1 Chromatogram of organic acids mixture solution

1: Phosphoric acid, 2: Citric acid, 3: Tartaric acid, 4: Malic acid, 5: Succinic acid, 6: Lactic acid, 7: Formic acid, 8: Acetic acid, 9: Pyroglutamic acid, 10: Propionic acid, 11: Isobutyric acid, 12: n-Butyric acid, 13: Isovaleric acid, 14: n-Valeric acid

Figure 2, 3, and 4 shows the chromatogram of rice wine, red wine, and plum wine. The red wine and plum wine are diluted ten times to measure the chromatograms.




Fig.2 Chromatogram of rice wine

1: Phosphoric acid, 2: Citric acid, 4: Malic acid, 5: Succinic acid, 6: Lactic acid, 8: Acetic acid, 9: Pyroglutamic acid




Fig.3 Chromatogram of red wine

1: Phosphoric acid, 2: Citric acid, 3: Tartaric acid, 4: Malic acid, 5: Succinic acid, 6: Lactic acid, 8: Acetic acid




Fig.4 Chromatogram of plum wine

1: Phosphoric acid, 2: Citric acid, 4: Malic acid, 5: Succinic acid, 8: Acetic acid