

## Analysis of Cyclodextrines by High Performance Liquid Chromatography with Evaporative Light Scattering Detection

### Introduction

Evaporative Light Scattering Detector (ELSD) is an universal detector for HPLC which detects the scattered light from the particles when illuminated of the involatile substance remained, after the sample is sprayed together with N<sub>2</sub> gas after the elution from column and heated to evaporate the mobile phase. LED is used as a light source to illuminate the involatile particles, and the light scattered will be collected and converted into the electronic signal by photomultiplier tube. Samples, such as sugar and fat being generally detected by UV absorption at short wavelength, or by using RI detector so far, can be measured with high sensitivity and stable baseline by using ELSD. Cyclodextrine is known as the oligosaccharide with its glucoses being connected in circle. Cyclodextrine is further named as  $\alpha$ -Cyclodextrine,  $\beta$ -Cyclodextrine and  $\gamma$ -Cyclodextrine for the cases that the number of glucoses becomes 6, 7, or 8. Since the property and the performance are different according to the number of glucoses, the amount of these cyclodextrines will influence the taste of food.

This report describes the analysis of cyclodextrine by using ELSD, polymer NH<sub>2</sub> column in HILIC mode.

**Keyword :** Cyclodextrine, HILIC, polymer NH<sub>2</sub> column, ELSD

### Experimental

#### Equipment


Pump: PU-2089  
 Autosampler: AS-2057  
 Column oven: CO-2060  
 Detector: ELS-2040

#### Conditions

Column: Shodex Asahipak NH2P-50 4E (4.6 mmID x 250 mmL)  
 Eluent: Water/Acetonitrile (40/60)  
 Flow rate: 1.0 mL/min  
 Column temp.: 50°C  
 ELSD condition: Nebulizer temp.: 30°C  
 Evaporator temp.: 30°C  
 Gas flow rate; 1.4 SLM  
 Injection volume: 10  $\mu$ L  
 Standard sample:  $\alpha$ ,  $\beta$ ,  $\gamma$ - Cyclodextrine 0.5 mg/mL each in Water/Acetonitrile (40/60)


### Result

Fig. 1 shows the chromatogram of  $\alpha$ ,  $\beta$ ,  $\gamma$ -Cyclodextrine, which were separated and detected properly.



**Fig. 1.** Chromatogram of the Standard mixture of Cyclodextrine  
 1:  $\alpha$ -Cyclodextrine , 2:  $\gamma$ -Cyclodextrine , 3:  $\beta$ -Cyclodextrine

Fig. 2 shows the chromatogram of green tea mixed with cyclodextrine. All the components were detected properly with the quantitative result :  $\alpha$ -Cyclodextrine 4.5 mg/10mL,  $\beta$ -Cyclodextrine 5.0 mg/10mL,  $\gamma$ -Cyclodextrine 22.7 mg/10mL.



**Fig. 2.** Chromatogram of green tea mixed with cyclodextrine

1:  $\alpha$ -Cyclodextrine , 2:  $\gamma$ -Cyclodextrine , 3:  $\beta$ -Cyclodextrine

Sample preparation: Green tea combined with cyclodextrine was mixed with mobile phase of the same amount, which was then filtrated by the membrane filter of 0.45 $\mu$ m.